İÇİNDEKİLER

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INDEX

PRODUCT USER MANUEL.....	33
PRODUCT DESCRIPTION.....	33
PURPOSE OF USE.....	33
TARGET PATIENT GROUP AND INTENDER USE.....	33
INDICATIONS.....	33
CONTRAINICATION.....	34
WARNINGS AND PRECAUTIONS.....	34
COMPLICATIONS AND SIDE EFFECTS.....	35
COMPETIBILITY INFORMATION.....	35
CLEANING AND DISINFECTING.....	35
STERILIZATION.....	35
PRODUCT MRI COMPETIBILITY INFORMATION.....	36
PROCEDURE.....	36
CLINICAL USAGE OF IMPLATECHONE ABUTMENT.....	36
PRODUCT USAGE.....	37
IMPLATECHONE BALL ABUTMENT.....	37
IMPLATECHONE LOCATOR ABUTMENT.....	39
IMPLATECHONE IMPLATECHONE 17 ° / 30 ° Multi Unit Abutment.....	41
IMPLATECHONE MULTI UNIT ABUTMENT.....	44
IMPLATECHONE CEMENT-RETAINED ABUTMENT.....	46
IMPLATECHONE TI-BASE ABUTMENT.....	50
IMPLATECHONE ABUTMENT AND PROSTHETIC SCREWS.....	53
IMPLATECHONE İYİLEŞME HEALING ABUTMENTS.....	57
DISPOSAL.....	60
DISCLAIMER OF LIABILITY.....	60
ADDITIONAL INFORMATION.....	60
SYMBOLS.....	61

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

## PRODUCT USER MANUAL

### 1. Product Description

The IMPLATECHONE prosthesis series consists of abutments in different types of diameters, lengths and platforms that are used in the restoration of IMPLATECHONE dental implants. They are available in a variety of shapes and sizes to suit the patient's needs.

IMPLATECHONE prostheses are used for the formation of connections between abutments and copings. IMPLATECHONE dental implants consist of different types, sizes and platforms.

### 2. Purpose of Use

IMPLATECHONE abutments are used to support dentures in order to fulfill the chewing function and to eliminate tooth deficiencies.

### 3. Target Patient Group And Intended User

IMPLATECHONE dental implants are intended to be used in patients with complete or partial edentulism who have completed their growth and development and do not have the conditions specified in the contraindications.

IMPLATECHONE dental implants may only be used by dentists.

Clinicians must have sufficient knowledge of implantology and practice skills, as well as comply with the instructions for use, in order to safely and properly use IMPLATECHONE dental implants.

### 4. Indications

They are screwed directly to the intraosseous implants to support the prosthesis.

- Complete and partially edentulous patients
- In chin and facial defects
- Complete tooth deficiencies characterized by excessively resorbed ridge
- In patients who have difficulties in the use of removable partial dentures
- In patients who refuse to use removable prosthesis
- In patients who want to have their teeth prepared
- In fixed prosthesis patients with very long spaces
- In severe damage to any edentulous area or soft tissues where the prosthesis fits
- Due to its high fracture strength, it can be used even in individuals with excessive overbite, bruxism or foreign body biting habits
- In cases where the height of the abutment is less than 7 mm and the axial thickness is less than 0.7 mm due to the patient's mouth closure, and in cases where the abutment must be angled more than 30 surgically and in the posterior region, it is higher to be used compared to ceramic abutments
- They can be used in the posterior shank and as a bridge abutment.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

## 5. Contraindication

Its use is contraindicated in individuals with hypersensitivity to titanium and titanium alloys. For implant treatment contraindications, refer to the relevant IMPLATECHONE implant instructions for use.

- Active Infection
- Foreign body sensitivity / material allergy (this situation should be determined and necessary precautions should be taken)
- Condition of bone structure and insufficient bone density (Performance of Dental Implant and Abutment Implant depends on bone condition and density)
- When the edentulous space is narrowed by the opposing and adjacent teeth
- In cases where the crown lengths of the abutment teeth are too short,
- In the presence of parafunctional habits such as bruxism
- When the use of a winged bridge (Kantilever) is designed
- Auxiliary teeth are not used if they lack adequate periodontal support.
- Inability of the corneal margin to adapt to the scallop structure of the gingiva;
- Deep gingival pockets formed as a result of not placing the implant deeper in order to provide a suitable aesthetic,
- Difficulty in cementation of the crown and removal of cement residues,
- The reflection in metallic blue color visible under the gingiva, aesthetic negativity in patients with thin gingival structure or high smile lines,
- Inadequate bone volume and/or quality
- Allergy or hypersensitivity to material

## 6. Warnings and Precautions

- Before the operation, the patient should be informed about the surgical risks and the positive and negative effects should be explained.
- The patient should be warned that the implanted device has a certain lifespan during which normal bone cannot replace, may be damaged by strenuous activity or trauma, and may be replaced in the future.
- It is recommended to review the product-specific surgical technique before the surgeon performs the surgical procedure. IMPLATECHONE dental implant can provide surgical technical information. Please contact the IMPLATECHONE dental implant sales representative.
- Choosing the right implant is extremely important. Appropriate type and size is taken into consideration, the patient's age, activity levels, bone density are preferred depending on whether or not he has undergone any previous surgical operation.
- IMPLATECHONE dental implant prostheses should be applied by securing respiration during oral use. Components or tools may injure the patient. Failure to follow these instructions may result in negative effect on breathing as swallowing the part or aspiration.
- IMPLATECHONE dental implant abutments should only be used with IMPLATECHONE dental implants on the appropriate platform.
- Abutments are disposable devices.
- Small diameter abutments should not be used where high mechanical loads are required.
- IMPLATECHONE dental implant prosthesis products should be used according to the instructions for use provided by the manufacturer. It is the physician's responsibility to decide whether the device is suitable for use according to these instructions and for the individual patient situation.
- For a successful implant treatment, it is important that the dentist and dental technician work in harmony.
- In order for the implant treatment to be successful in the long term, it is important that the patient's routine controls and oral hygiene training are given after the treatment is completed.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

## 7. Complications and Side Effects

The following reactions may occur during implant treatment:

- Hypersensitivity/allergic reactions
- Swallowing or aspiration of prosthetic components
- Peri-implantitis
- Implant component or prosthesis failure/fracture
- Poor aesthetic result
- Repetition of prosthetic treatment
- Minor bleeding
- Implant loss

## 8. Compatibility Information

All abutments can be used with the same abutment screwdriver.

There are different shapes of abutments in the IMPLATECHONE implant series. Abbreviations and color codes indicating the compatibility of implant and abutment platforms are available on the packages.

Implant Platform	Abbreviation	Color Code	Implant Diameter
Narrow Platform	NP	Violet	Ø 3.3 – Ø3.7
Regular Platform	RP	Turquoise	Ø4.1 – Ø4.7

## 9. Cleaning and Disinfecting

IMPLATECHONE abutments and parts are delivered unsterile. Before placing the products in the patient's mouth, they must be disassembled, cleaned and sterilized.

IMPLATECHONE recommends the following for cleaning and sterilization of abutments before use.

- Clean by brushing inside and outside with a brush under running water.
- The pre-treated product can be cleaned by hand, ultrasonic cleaner, or an automatic cleaning method.
- When selecting the automatic cleaning method, the appropriate detergent should be selected and the manufacturer's instructions should be followed.

## 10. Sterilization

IMPLATECHONE abutments and parts are delivered unsterile. IMPLATECHONE recommends the following for sterilization of abutments prior to use:

Method	Conditions	Drying
Humidity temperature (autoclave) pre-vacuum	121° C 30 min	Local

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

Caution: Use devices immediately after sterilization. Do not store sterilized devices.

## 11. Product MRI compatibility information

Non-clinical testing has demonstrated that Titanium Dental Materials are MR Compatible.

A patient with these system parts in his body can safely be scanned in MR systems if the following conditions are met:

Only 1.5 and 3.0 Tesla static magnetic fields

Magnetic spatial drop field up to 970 Gauss/cm ( 9.7 T/m) or less

Average specific absorption rate (SAR) for the whole body 2 W/kg (Normal Use Mode)

Under the above-mentioned scanning conditions, after 15 minutes of continuous scanning, our products are expected to generate a maximum temperature rise of 3°C.

## 12. Procedure

Be sure to maintain the prosthetic attachment by securing the abutment during polishing or other laboratory procedures.

- Insert the abutment into the implant analog on the working model.
- Make sure the analog-abutment joint holder elements are properly aligned.
- Fix the abutment by hand tightening the base screw.

### 12.1. Clinical Use of IMPLATECONE Abutment

The original abutment is taken from the laboratory. Remove the healing cap or temporary healer and rinse thoroughly then dry the inside of the implant. Remove the optimizer from the working model.

- Cleaning, disinfection and sterilization of dental healers described in chapter 9. And 10.
- The curative is placed in the mouth of patients.
- The screw is tightened after the abutment is properly placed.
- Be sure to fix the abutment implant with the appropriate screw.
- Tighten the abutment screw with the corresponding screwdriver.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

## 12.2. Product Usage

### 12.2.1. IMPLATECHONE Ball Abutment

IMPLATECHONE Ball abutments are manufactured using Ti6Al4V-ELI (ASTM F 136) material.

	Platform/H (mm)	NP-Narrow Platform		RP-Regular Platform	
		Ø3.3mm	Ø3.7mm	Ø4.1mm	Ø4.7mm
IMPLATECHONE Ball Abutment	H0,5	02.06.10.05		02.06.20.05	
	H1,0	02.06.10.10		02.06.20.10	
	H2,0	02.06.10.20		02.06.20.20	
	H3,0	02.06.10.30		02.06.20.30	
	H4,0	02.06.10.40		02.06.20.40	
	H6,0	02.06.10.60		02.06.20.60	

#### Intended Use

- IMPLATECHONE Ball abutments are used to support dentures in order to fulfill the chewing function and to eliminate tooth deficiencies. One-piece Ball abutments are screwed to Mode Medikal implants and form a basis for overdenture structures.

#### Indications

- They are screwed directly to the intraosseous implants to support the removable overdentures.

#### Contraindications

- Its use is contraindicated in individuals with hypersensitivity to titanium and titanium alloys. For implant treatment contraindications, refer to the relevant IMPLATECHONE implant instructions for use.

#### Complications And Side Effects

The following reactions may occur during implant treatment:

- Hypersensitivity/allergic reactions
- Swallowing or aspiration of prosthetic components
- Peri-implantitis
- Implant component or prosthesis failure/fracture
- Poor aesthetic result
- Repetition of prosthetic treatment
- Minor bleeding
- Implant loss

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

### Laboratory Procedure

**Note:** Screw the abutment to the analog during all laboratory procedures to preserve the abutment's attachment site to the implant.

- Combine the impression coping with the implant analogue.
- Obtain the plaster working model by using soft gingival silicone around the analogues.
- Determine the abutment gingival height according to the amount of gingiva in the relevant region in the plaster model obtained.
- Select the abutment suitable for the diameter of the dental implant used.
- Connect the selected abutment with suitable diameter and gingival height to the analog in the plaster model with the help of screws.
- Simply tighten the abutment screw by hand.
- Plan your prosthesis in accordance with the selected abutment type and fabricate it with the preferred fabrication method.
- Metal and plastic matrices used with ball abutments are compatible with Rhein 83.

### Clinical Procedure

- Remove the healing abutment located in the relevant area of the patient's mouth with the help of a screwdriver.
- Select the abutment at the appropriate gingival height, insert it into the related implant in the mouth and screw the abutment with 35Ncm torque. Torque values greater than 35 Ncm may damage the prosthetic parts. Torque values less than the recommended value may cause abutment loosening.
- If the prosthesis and metal housing are to be combined in a laboratory environment, take the impression with the relevant impression copings and custom impression tray.
- At the delivery of the prosthesis, select the appropriate retention plastic and place it in the metal housing.
- Place the prosthesis in the mouth and check the occlusal contacts, make the necessary adjustments.

**Oral and Prosthetic Care:** Good oral hygiene is extremely important for the success of the attachment. Ball attachments should be thoroughly cleaned daily to prevent plaque buildup, and the patient should use a soft, nylon-bristled toothbrush and non-abrasive toothpaste to clean the abutments. Plastic retention materials are subject to wear as part of normal use and may need to be replaced. Patients should be asked to continue routine follow-up visits for hygiene and attachment function evaluation. Follow-up visits are recommended at 6-month intervals. Patients should be examined for signs of inflammation around the implant abutments and for implant mobility.

**Placing and Removing Overdentures:** The patient should be given instructions on how to properly place the overdenture. The patient should ensure that the prosthesis is correctly placed on the abutments before applying pressure. Using both hands, the patient should press both sides of the overdenture prosthesis to firmly seat it.

The patient should not bite the overdenture prosthesis into place; otherwise the applied force will cause improper wear of the abutments, including the plastic retainers in the overdenture. To remove the overdenture, the patient should place their thumbs under the edges of the overdenture and remove the prosthesis by pushing from both sides.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

### 12.2.2. IMPLATECHONE Locator Abutment

IMPLATECHONE Locator abutments are manufactured using Ti6Al4V-ELI (ASTM F 136) material.

	Platform/H (mm)	NP - Narrow Platform		RP - Regular Platform	
		Ø3.3mm	Ø3.7mm	Ø4.1mm	Ø4.7mm
IMPLATECHONE Locator Abutment	H0,5	02.07.10.05		02.07.20.05	
	H1,0	02.07.10.10		02.07.20.10	
	H2,0	02.07.10.20		02.07.20.20	
	H3,0	02.07.10.30		02.07.20.30	
	H4,0	02.07.10.40		02.07.20.40	
	H6,0	02.07.10.60		02.07.20.60	

#### Intended use

- IMPLATECHONE Locator abutments are used to support dentures in order to fulfill the chewing function and to eliminate tooth deficiencies. One-piece Locator abutments are screwed to IMPLATECHONE implants and form a basis for overdenture structures.

#### Indications

- They are screwed directly to the intraosseous implants to support the removable overdentures.

#### Contraindications

- Its use is contraindicated in individuals with hypersensitivity to titanium and titanium alloys. For implant treatment contraindications, refer to the relevant IMPLATECHONE implant instructions for use.

#### Complications And Side Effects

The following reactions may occur during implant treatment:

- Hypersensitivity/allergic reactions
- Swallowing or aspiration of prosthetic components
- Peri-implantitis
- Implant component or prosthesis failure/fracture
- Poor aesthetic result
- Repetition of prosthetic treatment
- Minor bleeding
- Implant loss

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

### Laboratory Procedure

**Note:** Screw the abutment to the analog during all laboratory procedures to preserve the abutment's attachment site to the implant.

- Combine the impression coping with the implant analogue.
- Obtain the plaster working model by using soft gingival silicone around the analogues.
- Determine the abutment gingival height according to the amount of gingiva in the relevant region in the plaster model obtained.
- Select the abutment suitable for the diameter of the dental implant used.
- Connect the selected abutment with suitable diameter and gingival height to the analog in the plaster model with the help of screws.
- Simply tighten the abutment screw by hand.
- Plan your prosthesis in accordance with the selected abutment type and fabricate it with the preferred fabrication method.
- Metal and plastic matrices used with Locator abutments are compatible with Zest Locator and Kerator.

### Clinical Procedure

- Remove the healing abutment located in the relevant area of the patient's mouth with the help of a screwdriver.
- Select the abutment at the appropriate gingival height, insert it into the related implant in the mouth and screw the abutment with 35Ncm torque. Torque values greater than 35 Ncm may damage the prosthetic parts. Torque values less than the recommended value may cause abutment loosening.
- If the prosthesis and metal housing are to be combined in a laboratory environment, take the impression with the relevant impression copings and custom impression tray.
- At the delivery of the prosthesis, select the appropriate retention plastic and place it in the metal housing.
- Place the prosthesis in the mouth and check the occlusal contacts, make the necessary adjustments.

**Oral and Prosthetic Care:** Good oral hygiene is extremely important for the success of the attachment. Ball attachments should be thoroughly cleaned daily to prevent plaque buildup, and the patient should use a soft, nylon-bristled toothbrush and non-abrasive toothpaste to clean the abutments. Plastic retention materials are subject to wear as part of normal use and may need to be replaced. Patients should be asked to continue routine follow-up visits for hygiene and attachment function evaluation. Follow-up visits are recommended at 6-month intervals. Patients should be examined for signs of inflammation around the implant abutments and for implant mobility.

**Placing and Removing Overdentures:** The patient should be given instructions on how to properly place the overdenture. The patient should ensure that the prosthesis is correctly placed on the abutments before applying pressure. Using both hands, the patient should press both sides of the overdenture prosthesis to firmly seat it.

The patient should not bite the overdenture prosthesis into place; otherwise the applied force will cause improper wear of the abutments, including the plastic retainers in the overdenture. To remove the overdenture, the patient should place their thumbs under the edges of the overdenture and remove the prosthesis by pushing from both sides.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

### 12.2.3. IMPLATECHONE 17 ° / 30 ° Multi Unit Abutment

IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments are composed of two components. IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments are offered in two different angles, 17° and 30°.

IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments offers a wide range of prosthetic solutions such, bar-retained overdenture prostheses and screwed fixed restorations.

IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments covers and copings are manufactured using Ti6Al4V-ELI (ASTM F 136) material.

	Platform/H (mm)	NP - Narrow Platform		RP - Regular Platform	
		Ø3.3mm	Ø3.7mm	Ø4.1mm	Ø4.7mm
 17° Multi Base Abutment	H2,5	02.08.10.25		02.08.20.25	
	H3,5	02.08.10.35		02.08.20.35	
 30° Multi Base Abutment	H3,5	02.09.10.35		02.09.20.35	
	H4,0	02.09.10.40		02.09.20.40	
 Multi Base Ti Base Engaged Coping		03.02.05.01			
 Multi Base Ti Base Non Engaged Coping		03.02.05.02			

#### Intended use

- IMPLATECHONE abutments are used for support in eliminating tooth deficiencies in order to fulfill the chewing function. Two-piece IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments are attached to the implants with the help of screws and form the basis of prosthetic structures on implants such as screw retained fixed restorations.
- IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments Ti Base engaged and non engaged copings are supporting parts that provide a compatible interface with the IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments by cementing the screwed final restoration or temporary restoration in the laboratory environment.

#### Indications

- They are screwed directly to the intraosseous implants with the help of screws in order to support the implant-supported prostheses. They provide support for screw retained fixed restorations in edentulous or partially edentulous dental arches. IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments are used to connect bar-retained overdenture prostheses with the implant.

#### Contraindications

- Its use is contraindicated in individuals with hypersensitivity to titanium and titanium alloys. For implant treatment contraindications, refer to the relevant Mode Medical implant instructions for use.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

**Complications and side effects**

The following reactions may occur during implant treatment:

- Hypersensitivity/allergic reactions
- Swallowing or aspiration of prosthetic components
- Peri-implantitis
- Implant component or prosthesis failure/fracture
- Poor aesthetic result
- Repetition of prosthetic treatment
- Minor bleeding
- Implant loss

**Warning**

The abutment screw should be torqued with a torque value of 25 Ncm. Torque values greater than 25 Ncm may damage the prosthetic parts. Torque values less than the recommended value may cause abutment loosening. The prosthetic screws used in the IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments series should be torqued with a torque value of 15 Ncm.

**Clinical Procedure****Impression**

- IMPLATECHONE 17 ° / 30 ° Multi Unit Abutments are selected in accordance with implant location, angle and gingival height.
- The abutment body is placed in the appropriate position to the implant. The abutment carrier is flexible and can be bent if necessary.
- The abutment body is fixed to the implant with a torque ratchet. The abutment body should be torqued with a torque value of 25 Ncm.
- The abutment carrier is removed by turning.
- The abutment cover is placed on the abutment body and torqued with a standard hex wrench with a torque value of 25 Ncm.
- It is recommended to check the fit of the abutment by radiography.
- Implant impression is taken with open tray, closed tray or digital impression techniques.
- If provisional restoration will not be used, protective healing abutments are placed.
- If a provisional restoration is to be used, the provisional restoration is placed on the abutments and the prosthetic screw is tightened by hand.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

**Laboratory Procedure**

- Combine the impression coping with the implant analogue.
- Obtain the plaster working model by using soft gingival silicone around the analogues.
- Produce the prosthetic infrastructure according to the preferred method. Infrastructure can be produced as casting, CAD/CAM or metal sintering.
- Fabricate the final prosthesis according to the preferred laboratory technique.
- Metal housing and retention plastics used with Multi Base Locator Cover are compatible with Zest Locator and Kerator. Metal and plastic matrices used with Multi Base Ball Cover are compatible with Rhein 83.

**Screw Retained Fixed Restoration**

- Remove the healing abutment located in the relevant area of the patient's mouth in order to be able to try-in the prosthesis produced by the technician.
- If the patient has a provisional prosthesis with an occlusal screw, loosen the screws of the provisional prosthesis with the help of a screwdriver and remove the provisional prosthesis from the mouth.
- Loosen the prosthetic screw with a screwdriver to remove the prosthesis on the abutment in the plaster model.
- Remove the restoration from the model.
- Place the prosthesis on the relevant implants on the patient in the same position as the model and tighten the prosthetic screw with the help of a screwdriver and a torque ratchet. The prosthetic screw should be torqued with a torque value of 15 Ncm.

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#### 12.2.4. IMPLATECHONE Multi Unit Abutment

Multi Unit abutments are designed as one piece.

IMPLATECHONE Multi Unit abutments are manufactured using Ti6Al4V-ELI (ASTM F 136) material.

	Platform/H (mm)	NP - Narrow Platform		RP - Regular Platform	
		Ø3.3mm	Ø3.7mm	Ø4.1mm	Ø4.7mm
Multi-Unit Abutment	H0,5	02.05.10.05		02.05.20.05	
	H1,0	02.05.10.10		02.05.20.10	
	H2,0	02.05.10.20		02.05.20.20	
	H3,0	02.05.10.30		02.05.20.30	
	H4,0	02.05.10.40		02.05.20.40	

#### Intended Use

- IMPLATECHONE abutments are used for support in eliminating tooth deficiencies in order to fulfill chewing function. It is attached to the two-piece IMPLATECHONE implants with the help of screws and forms the basis for screw retained fixed implant restorations.

#### Indications

- They are attached directly to the intraosseous implants with the help of screws in order to support the implant-supported prostheses. They provide support for fixed prostheses in edentulous or partially edentulous dental arches.

#### Contraindications

- Its use is contraindicated in individuals with hypersensitivity to titanium and titanium alloys. For implant treatment contraindications, refer to the relevant IMPLATECHONE implant instructions for use.

#### Complications And Side Effects

The following reactions may occur during implant treatment:

- Hypersensitivity/allergic reactions
- Swallowing or aspiration of prosthetic components
- Peri-implantitis
- Implant component or prosthesis failure/fracture
- Poor aesthetic result
- Repetition of prosthetic treatment
- Minor bleeding
- Implant loss

#### Warning

- The Multi Unit abutment should be torqued to 25 Ncm. Torque values greater than 25 Ncm may damage the prosthetic parts. Torque values less than the recommended value may cause abutment loosening. The prosthetic screws used in the Multi Unit abutment series should be torqued with a torque value of 15 Ncm.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

## Clinical Procedure

### Impression

- The Multi Unit abutment is selected in accordance with the gingival height.
- The abutment is attached and tightened to the implant with a torque ratchet. The abutment should be torqued with a torque value of 25 Ncm.
- It is recommended to check the fit of the abutment by radiography.
- Implant impression is taken with open tray, closed tray or digital impression techniques.
- If provisional restoration will not be used, protective healing abutments are placed.
- If a provisional restoration is to be used, the provisional restoration is placed on the abutments and the prosthetic screw is tightened by hand.

### Laboratory Procedure

- Combine the impression transfer with the implant analogue.
- Obtain the plaster working model by using soft gingival silicone around the analogues.
- Produce the prosthetic infrastructure according to the preferred method. Infrastructure can be produced as casting, CAD/CAM or metal sintering.
- Fabricate the final prosthesis according to the preferred laboratory technique.

### Restoration Delivery

- Remove the healing abutment located in the relevant area of the patient's mouth in order to be able to try-in the prosthesis produced by the technician.
- If the patient has a provisional prosthesis with an occlusal screw, loosen the screws of the provisional prosthesis with the help of a screwdriver and remove the provisional prosthesis from the mouth.
- Loosen the prosthetic screw with a screwdriver to remove the prosthesis on the abutment in the plaster model.
- Remove the restoration from the model.
- Place the prosthesis on the relevant implants on the patient in the same position as the model and tighten the prosthetic screw with the help of a screwdriver and a torque ratchet. The prosthetic screw should be torqued with a torque value of 15 Ncm.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

### 12.2.5. IMPLATECHONE Cement-Retained Abutments

IMPLATECHONE cement-retained abutments are manufactured using Ti6Al4V-ELI (ASTM F 136) material.

	Platform/H (mm)	NP - Narrow Platform		RP - Regular Platform	
		Ø3.3mm	Ø3.7mm	Ø4.1mm	Ø4.7mm
 Profil Abutment	H1,0	02.10.10.10		02.10.20.10	
	H2,0	02.10.10.20		02.10.20.20	
	H3,0	02.10.10.30		02.10.20.30	
 Aesthetic Abutment	H1,0	02.04.10.10		02.04.20.10	
	H2,0	02.04.10.20		02.04.20.20	
	H3,0	02.04.10.30		02.04.20.30	
 15° Angled Abutment	H1,0	02.02.10.10		02.02.20.10	
	H2,0	02.02.10.20		02.02.20.20	
	H3,0	02.02.10.30		02.02.20.30	
	H4,0	02.02.10.40		02.02.20.40	
 25° Angled Abutment	H1,0	02.03.10.10		02.03.20.10	
	H2,0	02.03.10.20		02.03.20.20	
	H3,0	02.03.10.30		02.03.20.30	
	H4,0	02.03.10.40		02.03.20.40	
 Premill Abutment		03.01.10.01		03.01.20.01	

#### Intended Use

- IMPLATECHONE abutments are used to support fixed implant supported prostheses in order to fulfill the chewing function and to eliminate tooth deficiencies. It is attached to the two-piece IMPLATECHONE implants with the help

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of abutment screws and forms a basis for cement retained prosthetic structures such as crowns and bridges.

		NP – Narrow Platform		RP - Regular Platform	
Platform/H (mm)		Ø3.3mm	Ø3.7mm	Ø4.1mm	Ø4.7mm
Abutment çapı		Ø3,5mm	Ø4,0mm	Ø4,5mm	Ø5,0mm
 Dual Abutment	H0,5	02.01.35.05	02.01.40.05	02.01.45.05	02.01.50.05
	H1,0	02.01.35.10	02.01.40.10	02.01.45.10	02.01.50.10
	H2,0	02.01.35.20	02.01.40.20	02.01.45.20	02.01.50.20
	H3,0	02.01.35.30	02.01.40.30	02.01.45.30	02.01.50.30
	H4,0	02.01.35.40	02.01.40.40	02.01.45.40	02.01.50.40

### Indications

- They are attached directly to the intraosseous implants with the help of abutment screws in order to support the implant-supported cement retained fixed prostheses.

Abutment type	Retention	Indications	
		Crown	Bridge
Profile Abutment	Cemented	✓	✓
Dual Abutment	Cemented	✓	✓
Aesthetic Abutment	Cemented	✓	✓
15°/25° Angled Abutment	Cemented	✓	✓
Premill Abutment	Cemented	✓	✓

### Contraindications

- Its use is contraindicated in individuals with hypersensitivity to titanium and titanium alloys. For implant treatment contraindications, refer to the relevant Mode Medikal implant instructions for use.

### Complications And Side Effects

The following reactions may occur during implant treatment:

- Hypersensitivity/allergic reactions
- Swallowing or aspiration of prosthetic components
- Peri-implantitis
- Implant component or prosthesis failure/fracture
- Poor aesthetic result
- Repetition of prosthetic treatment
- Minor bleeding
- Implant loss

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**Warning**

- Before the delivery of the prosthesis, the abutment placed in the patient's mouth should be torqued with a torque value of 35 Ncm. Torque values greater than 35 Ncm may damage the prosthetic parts. Torque values less than the recommended value may cause abutment loosening.

**Laboratory Procedure****Traditional Workflow in Implant Supported Fixed Dentures**

**Note:** Screw the abutment to the analog during all laboratory procedures to preserve the abutment's attachment to the implant.

- Combine the impression coping with the implant analogue.
- Obtain the plaster working model by using soft gingival silicone around the analogues.
- Determine the abutment gingival height according to the amount of gingiva in the relevant region in the plaster model obtained.
- Select the abutment suitable for the diameter of the dental implant used.
- Connect the selected abutment with suitable diameter and gingival height to the analog in the plaster model with the help of screws.
- Simply tighten the abutment screw by hand.
- Plan your prosthesis in accordance with the selected abutment type and fabricate it with the preferred fabrication method.

**Digital Workflow in Implant Supported Fixed Dentures**

- Digital prosthetic workflow can be carried out by using STL data obtained with intraoral scanners or by obtaining STL data of impressions taken with traditional methods with a desktop scanner.
- The use of software that includes IMPLATECHONE implants in their library for prosthesis design and production on the digital platform will ensure that the internal geometry of the abutment and restoration connection is matching (eg Exocad DentalCAD).

**Clinical Procedure for Implant Supported Fixed Prostheses**

- For the try-in of the prosthesis produced by the technician, remove the healing abutment located in the relevant area of the patient's mouth with the help of a screwdriver.
- Loosen the abutment screw with a screwdriver to remove from the model.
- Remove the abutment from the model.
- Place the abutment on the patient's implant in the same position as the model and tighten the screw with the help of a screwdriver.
- Perform the prosthesis try-in.
- After try-in, remove the prosthesis from the patient's mouth and loosen the abutment screw with the help of a screwdriver.
- Place the abutment on the corresponding analog on the model and tighten the abutment screw with the help of the screwdriver.
- Before cementation of the prosthesis, tighten the screw of the abutment placed in the relevant area of the patient's mouth with a torque value of 35 Ncm.
- It is recommended to check the fit of the abutment by radiography.
- Cover the screw hole with a suitable material and cement the prosthesis with the preferred luting cement.

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**Warning:** After the abutment screw is torqued, the abutment screw hole should be closed with suitable materials. Correct closure of the screw hole makes it possible to access the screw in cases where the abutment needs to be separated from the implant.

#### Abutment Modification

When necessary, abutments can be modified according to the patient's anatomy by the technician in the dental laboratory or by the dentist in the clinic.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

### 12.2.6. IMPLATECHONE Ti-Base Abutments

IMPLATECHONE Ti-base abutments are manufactured using Ti6Al4V-ELI (ASTM F 136) material.

	Platform/H (mm)	NP - Narrow Platform		RP - Regular Platform	
		Ø3.3mm	Ø3.7mm	Ø4.1mm	Ø4.7mm
Ti-Base Engaged (Digital) Abutment	H0,7	03.02.10.07		03.02.20.07	
	H2,5	03.02.10.25		03.02.20.25	
Ti-Base Non-Engaged (Digital) Abutment	H0,7	03.01.10.07		03.01.20.07	
	H2,5	03.01.10.25		03.01.20.25	

#### Intended use

- IMPLATECHONE Ti-base abutments are used to support fixed implant supported prostheses in order to fulfill the chewing function and to eliminate tooth deficiencies. IMPLATECHONE implants with the help of abutment screws and forms a basis for screw retained prosthetic structures such as crowns and bridges.

#### Indication

- They are attached directly to the intraosseous implants with the help of abutment screws in order to support the implant-supported cement or screw retained fixed prostheses. IMPLATECHONE Ti-base abutments are titanium bases that can also be used in personal abutment construction.

Abutment type	Retention	Indications	
		Crown	Bridge
Ti-Base Engaged (Digital) Abutment	Screw	✓	
Ti-Base Non-Engaged (Digital) Abutment	Screw		✓

#### Contraindications

- Its use is contraindicated in individuals with hypersensitivity to titanium and titanium alloys. For implant treatment contraindications, refer to the relevant Mode Medikal implant instructions for use.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

**Complications And Side Effects**

The following reactions may occur during implant treatment:

- Hypersensitivity/allergic reactions
- Swallowing or aspiration of prosthetic components
- Peri-implantitis
- Implant component or prosthesis failure/fracture
- Poor aesthetic result
- Repetition of prosthetic treatment
- Minor bleeding
- Implant loss

**Warning**

- Before the delivery of the prosthesis, the abutment placed in the patient's mouth should be torqued with a torque value of 35 Ncm. Torque values greater than 35 Ncm may damage the prosthetic parts. Torque values less than the recommended value may cause abutment loosening.

**Procedure****Laboratory Procedure**

- Combine the impression coping with the implant analogue.
- Obtain the plaster working model by using soft gingival silicone around the analogues.
- Determine the abutment gingival height according to the amount of gingiva in the relevant region in the plaster model obtained.
- Select the abutment suitable for the diameter of the dental implant used.
- Connect the selected abutment with suitable diameter and gingival height to the analog in the plaster model with the help of screws.
- Simply tighten the abutment screw by hand.
- Plan your prosthesis with screw or cement retained in accordance with the selected abutment type and produce it with the preferred production method.
- Place the fabricated restoration to the analog before cementing it to the Ti-base abutment.
- Seal the screw channel with wax.
- Apply a self-curing resin cement to the abutment. Use the resin cement in accordance with the manufacturer's instructions.
- Place the restoration on the Ti-base abutment.
- Immediately clean any cement overflowing from the abutment. Polish the joint after the cement has hardened.

**Digital Workflow in Implant Supported Fixed Dentures**

- Digital prosthetic workflow can be carried out by using STL data obtained with intraoral scanners or by obtaining STL data of impressions taken with traditional methods with a desktop scanner.
- The use of software that includes IMPLATECHONE implants in their library for prosthesis design and production on the digital platform will ensure that the internal geometry of the abutment and restoration connection is matching (eg Exocad DentalCAD).

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

**Clinical Procedure for Implant Supported Fixed Protheses**

- For the try-in of the prosthesis produced by the technician, remove the healing abutment located in the relevant area of the patient's mouth with the help of a screwdriver.
- Loosen the abutment screw with a screwdriver to remove from the model.
- Remove the abutment from the model.
- Place the abutment on the patient's implant in the same position as the model and tighten the screw with the help of a screwdriver.
- Perform the prosthesis try-in.
- After try-in, remove the prosthesis from the patient's mouth and loosen the abutment screw with the help of a screwdriver.
- Place the abutment on the corresponding analog on the model and tighten the abutment screw with the help of the screwdriver.
- Before cementation of the prosthesis, tighten the screw of the abutment placed in the relevant area of the patient's mouth with a torque value of 35 Ncm.
- It is recommended to check the fit of the abutment by radiography.
- Cover the screw hole with a suitable material and cement the prosthesis with the preferred luting cement.

**Warning:** After the abutment screw is torqued, the abutment screw hole should be closed with suitable materials. Correct closure of the screw hole makes it possible to access the screw in cases where the abutment needs to be separated from the implant.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

### 12.2.7. IMPLATECHONE Abutment and Prosthetic Screws

IMPLATECHONE abutment and prosthetic screws are part of IMPLATECHONE Dental Implant Systems, which is an integrated system with IMPLATECHONE dental implants and related abutments, cover screws, healing abutments, surgical and prosthetic parts and instruments; It is produced as a single piece using Ti6Al4V-ELI (ASTM F 136). IMPLATECHONE abutment screws are implant screws that allow abutments to be fixed to implants. IMPLATECHONE prosthetic screws are implant screws that enable the fixation of patient-specific temporary or permanent crowns to the relevant abutments.

The abutments that IMPLATECHONE abutments and prosthetic screws are compatible with are shown in Table 1 below.

**Table 1**

Screw	Abutment	Screwdriver
Abutment Screws	Dual Abutment	IMPLATECHONE Hex Screwdriver
	Aesthetic Abutment	
	Angled Abutment 15°	
	Angled Abutment 25°	
	Profile Abutment	
	Ti-base Non-Engaged Abutment	
	Ti-base Engaged Abutment	
	Premill Abutment	
Prosthetic Screws	Multi Unit Abutment	IMPLATECHONE Hex Screwdriver
	Multi Base Abutment	

PROSTHETIC SCREW	
Platform	NP/RP
REF NO	03.02.05.03

ABUTMENT SCREW	
Platform	NP/RP
REF NO	03.02.15.01

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

**Intended Use**

- IMPLATECHONE abutment screws are intended to be used to fix IMPLATECHONE abutments to corresponding IMPLATECHONE dental implants. IMPLATECHONE prosthetic screws are intended to be used to fix patient-specific temporary or permanent crowns to the relevant abutments.

**Indications**

- IMPLATECHONE abutment screws are indicated to be used to fix the abutments which are related to the prostheses to be produced for the functional, phonetic and aesthetic rehabilitation of the upper and/or lower jaw in totally or partially edentulous patients, to the placed IMPLATECHONE dental implants. IMPLATECHONE prosthetic screws are indicated for use to fix prostheses to IMPLATECHONE abutments designed for the functional, phonetic and aesthetic rehabilitation of the upper and/or lower jaw in totally or partially edentulous patients.

**Contraindications**

The use of IMPLATECHONE abutment and prosthetic screws is contraindicated in the presence of the following conditions;

- Patients medically unfit for implant treatment
- Allergy or hypersensitivity to Ti6Al4V-ELI (ASTM F 136) materials

Note: For contraindications to implant treatment, refer to the relevant IMPLATECHONE implant instructions for use; for the contraindications of the abutments where the screws will be used, refer to the relevant IMPLATECHONE abutment instructions for use.

**Complications And Side Effects**

Dental implant treatment outcomes can be affected by many variables. Possible complications and side effects listed below may be observed after the use of IMPLATECHONE abutment and prosthetic screws.

- Local pain
- Micro hemorrhages
- Swelling
- Local inflammations
- Gingival injuries
- Hypersensitivity and/or allergic reactions
- Other toxicity reactions
- Aspiration or swallowing of parts during operation
- Triggering of the pharyngeal reflex during the insertion of the product
- Screw fracture
- Screw loosening
- Component or prosthesis failure/fracture

**Warnings**

- IMPLATECHONE abutment and prosthetic screws are intended to be use with IMPLATECHONE dental implants and related abutments in patients with complete or partial edentulism who have completed growth and development and do not have the conditions specified in the contraindications.
- IMPLATECHONE abutment and prosthetic screws should only be used in prosthesis stages with their own system elements and compatible abutments. The use of products with different brands and materials may lead to mechanical problems, failure of implants, tissue damage or aesthetic dissatisfaction.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

- IMPLATECHONE abutment and prosthetic screws are packaged non-sterile. It is the clinician's responsibility to ensure sterility prior to use.
- Sterile handling is essential. Never use potentially contaminated components. Contamination can cause infections.
- IMPLATECHONE abutments and prosthetic screws are for single use and should not be reused.
- The recommended torque value is 25 Ncm for abutment screws and 15 Ncm for prosthetic screws.
- All surgical instruments and tooling used must be kept in good condition and care must be taken not to damage the implants or other components of the instruments.

**Pre-operative**

- A comprehensive clinical and radiological examination is necessary to determine the psychological and physical condition of the patient.
- Preoperative hard and soft tissue deficiencies may result in undesirable aesthetic results and/or negative implant angulation.
- Particular attention should be paid to patients with local or systemic factors that may affect hard tissue healing, soft tissue healing, or osseointegration of implanted implants. (eg Type 1 Diabetes, bone metabolism diseases, bleeding disorders, anticoagulant therapy, parafunctional habits such as bruxism, smoking, poor oral hygiene, pregnancy, acute periodontal diseases, orofacial radiotherapy, adjacent tissue infections, etc.)
- Particular attention should be paid to patients receiving bisphosphonate treatment.
- Implant treatment should not be applied to patients who have not completed bone growth and development.
- Placement of implants and design of prostheses should be specific to the patient. In the presence of extreme bruxism and other parafunctional habits and jaw relationship disorders, different treatment options can be planned.

**During the operation**

- Before placing the IMPLATECHONE abutment and prosthetic screws, ensure that the inner surface of the implant is clean and free from blood.
- Be careful against the risk of aspiration and swallowing due to the size of the products during intraoral application. Aspiration of products can lead to infection or unwanted physical injury.

**Post-operative**

- For the long-term success of treatments with the IMPLATECHONE dental implant system, it is important to follow up the patient regularly and to give the necessary oral hygiene training to the patient.

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### Procedure

- Use the implant screw suitable for the abutment you selected for treatment. (See Table 1)
- Insert the IMPLATECHONE implant screw into the IMPLATECHONE screwdriver.
- First, tighten the screw you placed on the abutment by hand.
- The IMPLATECHONE torque ratchet and IMPLATECHONE screwdriver should be used to fix the implant screws.

**CAUTION:** When using IMPLATECHONE abutment and prosthetic screws, it is recommended to ensure that the screwdriver is firmly seated on the screws to avoid aspiration/swallowing of the screws.

**NOTE:** The recommended torque value for tightening is 25 Ncm for abutment screws and 15 Ncm for prosthetic screws.

**CAUTION:** Torque of IMPLATECHONE abutment and prosthetic screws higher than the recommended torque values may cause screw fractures, and under-torque may cause loosening of the screws.

### Compatibility Information

IMPLATECHONE abutment and prosthetic screws are compatible with IMPLATECHONE implant systems products and components. Do not use implants and components made of different metals and supplied from different manufacturers with IMPLATECHONE dental implant system products. IMPLATECHONE abutments and prosthetic screws can only be used with the corresponding abutments.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

## 12.2.8. IMPLATECHONE HEALING ABUTMENTS

**NP**  
**Healing**  
**Abutments**

**RP**  
**Healing**  
**Abutments**

IMPLATECHONE healing abutments are part of IMPLATECHONE Dental Implant Systems, which is an integrated system with IMPLATECHONE dental implants and related abutments, cover screws, surgical and prosthetic parts and instruments; It is produced as a single piece using Ti6Al4V-ELI (ASTM F 136) material.

- IMPLATECHONE healing abutments are components that support gingival healing on the implant, provide soft tissue contouring, and are directly connected to endosseous dental implants.
- IMPLATECHONE healing abutments are designed to be compatible with all the following IMPLATECHONE dental implants in three different platforms, narrow platform (NP) and regular platform (RP)

Implant / Platform	Narrow Platform (NP)	Regular Platform (RP)
IMPLATECHONE Implant	NP Healing Abutments	RP Healing Abutments
IMPLATECH ONE-S SHORT Implant		

Diameter and length options of IMPLATECHONE healing abutments are shown in the table below.

Platform	NP	NP	RP	RP
Length (L)	Ø4	Ø4.5	Ø4.5	Ø5
2 mm	03.40.10.02	03.45.10.02	03.45.20.02	03.50.20.02
4 mm	03.40.10.04	03.45.10.04	03.45.20.04	03.50.20.04
6 mm	03.40.10.06	03.45.10.06	03.45.20.06	03.50.20.06

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

**Intended Use**

- They are intended to be used for supporting gingival healing on the implant and providing soft tissue contouring by placing them on the corresponding IMPLATECHONE dental implants.

**Indications**

- IMPLATECHONE healing abutments are indicated for use with IMPLATECHONE dental implants placed for the functional, phonetic and aesthetic rehabilitation of the upper and/or lower jaw in totally or partially edentulous patients.

**Contraindications**

The use of IMPLATECHONE healing abutments is contraindicated in the presence of the following conditions;

- Patients who are medically unfit for implant treatment
- Allergy or hypersensitivity to Ti6Al4V-ELI (ASTM F 136) materials

**Complications And Side Effects**

Dental implant treatment outcomes can be affected by many variables. Possible complications and side effects listed below may be observed after the use of Mode healing abutments.

- Local pain
- Micro hemorrhages
- Swelling
- Local inflammations
- Gingival injuries
- Hypersensitivity and/or allergic reactions
- Other toxicity reactions
- Aspiration or swallowing of parts during operation
- Triggering of the pharyngeal reflex during the insertion of the product
- Growth of soft tissue on the healing abutments during the healing process
- Calculus formation on the healing abutments

**Warning**

- IMPLATECHONE healing abutments are intended for use with IMPLATECHONE dental implants in patients with complete or partial edentulism who have completed growth and development and do not have the conditions specified in the contraindications.
- IMPLATECHONE healing abutments should only be used with their own system elements and paired with the correct platform during the surgical and prosthetic phases. The use of products with different brands and materials may lead to mechanical problems, failure of implants, tissue damage or aesthetic dissatisfaction.
- IMPLATECHONE healing abutments should be used only during the healing process. They cannot be used to support a restoration.
- IMPLATECHONE healing abutments are packaged non-sterile. It is the clinician's responsibility to ensure sterility prior to use.
- Sterile handling is essential. Never use potentially contaminated components. Contamination can cause infections.
- All surgical instruments and tooling used must be kept in good condition and care must be taken not to damage the implants or other components of the instruments.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

**Pre-operative**

- A comprehensive clinical and radiological examination is necessary to determine the psychological and physical condition of the patient.
- Preoperative hard and soft tissue deficiencies may result in undesirable aesthetic results and/or negative implant angulation.
- Particular attention should be paid to patients with local or systemic factors that may affect hard tissue healing, soft tissue healing, or osseointegration of implanted implants. (eg Type 1 Diabetes, bone metabolism diseases, bleeding disorders, anticoagulant therapy, parafunctional habits such as bruxism, smoking, poor oral hygiene, pregnancy, acute periodontal diseases, orofacial radiotherapy, adjacent tissue infections, etc.)
- Particular attention should be paid to patients receiving bisphosphonate treatment.
- Implant treatment should not be applied to patients who have not completed bone growth and development.
- Placement of implants and design of prostheses should be specific to the patient. In the presence of extreme bruxism and other parafunctional habits and jaw relationship disorders, different treatment options can be planned.

**During the operation**

- Before placing the IMPLATECHONE healing abutments, ensure that the inner surface of the implant is clean and free from blood.
- Be careful against the risk of aspiration and swallowing due to the size of the products during intraoral application. Aspiration of products can lead to infection or unwanted physical injury.

**Post-operative**

- For the long-term success of treatments with the IMPLATECHONE dental implant system, it is important to follow up the patient regularly and to give the necessary oral hygiene training to the patient.

**Procedure**

- Select the appropriate healing abutment for the implant platform and gingival height.
- Take the IMPLATECHONE healing abutment with the IMPLATECHONE screwdriver.
- Connect the healing cap to the implant you placed and tighten it by hand.
- To remove the IMPLATECHONE healing abutment, loosen it by hand using the screwdriver.

**CAUTION:** When using IMPLATECHONE healing abutments, it is recommended to ensure that the screwdriver is firmly seated in the head to avoid the risk of aspiration/swallowing of the head.

**NOTE:** The recommended tightening torque is 5-10 Ncm manually by hand.

**Compatibility Information**

IMPLATECHONE healing abutments are compatible with IMPLATECHONE implant systems products and components. Do not use implants and components made of different metals and supplied from different manufacturers with IMPLATECHONE dental implant system products. IMPLATECHONE healing abutments can only be used with those with the corresponding connection on all IMPLATECHONE dental implants.

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR

### 13. Disposal

Disposal should be done in an environmentally sustainable manner according to local regulations. Hazardous waste and cutting/drilling tools for contaminated devices should be disposed of in suitable containers that meet special technical requirements.

### 14. Disclaimer of Liability

These products are part of a general concept and should only be used with the relevant original products in accordance with IMPLATECHONE's instructions and recommendations. The use of products not distributed by IMPLATECHONE and made by third parties will void any expressed or implied warranty or other obligations of IMPLATECHONE.

IMPLATECHONE products should be used in accordance with the instructions for use provided by the manufacturer. It is the clinician's responsibility to use products in accordance with these instructions and to determine whether the product is appropriate for the individual patient situation. The clinician is also responsible for regularly following the latest updates on the relevant IMPLATECHONE product and its applications. IMPLATECHONE is not responsible for any direct, indirect, penal or other damages arising from or in connection with any errors in professional judgment or practice related to the use of IMPLATECHONE products; assumes no liability, express or implied.

### 15. Additional Information

Non Sterile products do not have a specified shelf life.

#### Life Time

The useful life of the products is determined as 10 years.

#### NOTICE REGARDING SERIOUS INCIDENT

For a patient/user/third party in the European Union and in countries with an identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. The contact information for the manufacturer of this device to report a serious incident is as follows:

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Onaylayan  
ŞÜKRİYE AKARLAR

## 16. Symbols

ISO 15223-1: 2012 Symbols and Meanings Prepared According to the Standard

2292

Notified Body Number

Non Sterile

Attention

Lot number

Do Not Use For The Second  
Time

Barcode

See User Manual

Manufacturer Information


Reference No

## Product Shelf Life Information

The product does not have a specified shelf life. The service life of the products has been determined as 10 years with the support of the literature.

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Symbol	Symbol Description
	Patient name or patient ID
	Date of implantation
	Name and Address of the implanting healthcare institution/provider
	Name and Address of the manufacturer
	Information website for patients
	Device name
	Lot Number/Batch Code
	Explanation of unique device identifier as AIDC Format <i>AIDC: Automatic identification and data capture format (e.g. linear or 2D-Barcodes)</i>

Aksaray Mah. Hobyar Mektebi Sok. No.49/1 Fatih – İstanbul, Türkiye

Hazırlayan  
ESEN ERDOĞAN

Onaylayan  
ŞÜKRİYE AKARLAR