**KULLANMA KILAVUZU /
USER MANUAL**

İÇİNDEKİLER

IMPLATECH ONE IMPLANT.....	2
IMPLATECH ONE-S SHORT IMPLANT.....	8
IMPLATECH ONE İMPLANT KAPATMA VİDALARI.....	14

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approverd By
ŞÜKRİYE AKARLAR

**KULLANMA KILAVUZU /
USER MANUAL**

IMPLATECH ONE İMPLANT

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IMPLATECH ONE IMPLANT

1. Product Description

Implatech One Dental Implants are manufactured from biocompatible titanium or titanium alloy. Implatech One Dental Implants contain various surface treatments. Please refer to individual product labels for specific product description.

Implatech One dental implants are part of Implatech One Dental Implant Systems, which is an integrated system with related abutments, healing abutments, cover screws, surgical and prosthetic parts and instruments; They are endosseous screw type implants that have an internal conical octagon connection between the implant and the abutment.

Surface roughening of ImplatechOne implants produced from biocompatible Titanium Grade 4 (ASTM F 67) material is done with biphasic calcium phosphate (BCP).

They are packaged together with the cover screw produced using Ti6Al4V-ELI (ASTM F 136) material.

Implatech One implants are produced in three different platforms as narrow platform (NP), regular platform (RP)

Diameter and length options of Implatech One implants according to their platforms are given in the table below.

Platform	NP	NP	RP	RP
Implant Ø(D)	Ø3.3	Ø3.7	Ø4.1	Ø4.7
Uzunluk (L)				
8 mm	01.01.33.08	01.01.37.08	01.01.41.08	01.01.47.08
10 mm	01.01.33.10	01.01.37.10	01.01.41.10	01.01.47.10
11,5 mm	01.01.33.11	01.01.37.11	01.01.41.11	01.01.47.11
13 mm	01.01.33.13	01.01.37.13	01.01.41.13	01.01.47.13

2. Intended Use

IMPLATECH ONE IMPLANT dental implants;

The implant, which is the main part of the dental implant system, is known as the artificial tooth root. It adheres effectively to the jawbone with its sandblasted and hydrophilic outer surface. Parts of Dental Implant Systems can be used to support and/or retain single-member, multi-member or all-arch fixed and/or removable dentures in partially or fully edentulous maxillary or mandibular arches. Implatech One implants are intended to be used in oral implantation as a support for prosthetic applications to be made in order to eliminate the function, phonation and aesthetic problems caused by tooth deficiencies in the lower and/or upper jaw bones.

3. Target Patient Group And Intended User

Implatech One implants are intended to be used in patients with complete or partial edentulism who have completed their growth and development and do not have the conditions specified in the contraindications.

Implatech One dental implants are intended to be used by dentists only.

Clinicians must have sufficient knowledge of implantology and practice skills, as well as comply with the instructions for use, in order to use Implatech One dental implants safely and properly.

4. Indications

- Complete and partially edentulous patients
- Jaw and facial defects

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

**KULLANMA KILAVUZU /
USER MANUAL**

- In complete tooth deficiencies with excessively resorbed ridges
- In patients who have difficulties in the use of removable partial dentures
- In fixed prosthesis patients with very long spaces
- In patients who refuse to use removable prosthesis
- In patients who do not want to have their teeth prepared
- In cases of serious changes in any edentulous area or in the soft tissues where the full denture sits
- In cases where oral muscular coordination is impaired
- In cases where tissue tolerance is low
- In patients with parafunctional habits that destabilize the prosthesis
- In patients with high expectations from full dentures
- In patients with excessive vomiting reflex
- In patients who are psychologically opposed to the removable prosthesis
- In the presence of insufficient number and position of supporting teeth
- In the absence of a single tooth with healthy adjacent teeth
- Trauma-related tooth loss or root fractures
- In the presence of internal granulomas
- In the presence of non-vital teeth with apical periodontitis that cannot be corrected by conservative methods or surgery
- In external and internal root resorptions
- For orthodontic anchorage purposes
- Dental agenesis
- Conservative treatment request (the patient's desire not to interfere with his healthy teeth)

5. Contraindication

- Recent myocardial infarction. The first thing to consider here is surgical stress. Because the adrenaline that will be released during a surgical procedure will adversely affect the cardiac condition of the patient. In addition, there may be uncontrollable vasoconstriction, which may cause some disturbances in heart rhythm. Since almost all of these errors use anticoagulant drugs, a coagulation disorder may be in question. There is also a risk of infections.
- Those with a heart valve prosthesis. Likewise, there is a risk against surgical stress, coagulation imbalance and infections leading to the loss of the heart valve.
- Severe kidney diseases. Here, calcium is not absorbed from the tubules. This causes metabolic loss of calcium. Parathormone dysfunction. It causes metabolic osteopenia. risk of infection.
- Treatment of severe osteomalacia (rickets) There is hypophosphocalcic bone and osteoidosis. Failure to integrate the implant with the bone. 75% of rickets receive vit D + Ca. This hinders integration. There is a risk of infection.
- Generalized secondary osteoporosis. Although primary osteoporosis is physiological, secondary generalized osteoporosis has pathological conditions such as Hodking's disease. Rarefaction is seen in the bone structure. Absence of osteodosis; Although there is a mass increase in the bone, there is a volumetric decrease. Integration of the bone with the implant does not occur. There is a risk of infection.
- Uncontrolled Diabetes Mellitus. Hyperosmolarity is seen in the blood. Dehydration and metabolic diseases are in question. Angiopathy: Most diabetics suffer from micro and macro angiopathy. This leads to a predisposition to tissue degeneration. There is a risk of infection. In addition, wound healing is difficult.
- Those undergoing radiotherapy. The defense mechanism is broken. Osteoinduction and osteoconduction are impaired. Physiological periosteal activity is impaired. There is a predisposition to tissue necrosis and a risk of infection.
- Chronic or severe alcoholism. Liver diseases such as cirrhosis are in question. These, in turn, eventually lead to coagulation disorders. Medullary disorders (vit B1, B6, B12). As a result of these, there may be risk of anemia, thrombocyte disorders, hemorrhage in different distant places. Delay in recovery occurs due to malnutrition. Psychological disorders. There is a risk of infection.
- Severe hormonal disorders. There is a metabolic calcium deficiency. There is deterioration of the implant bed.
- Drug addiction. Loss of sense of pruritis Nutritional disorder. Loss of resistance to diseases. Psychological disorders. There is a risk of infection.
- Long-term use of immunosuppressant drugs. Delay in recovery. Medullary aplasia. Increased bone fragility. risk of infection
- Severe connective tissue diseases. For example, those that cause permanent damage to connective tissue, such as Lupus Erythematosus.
- Severe blood diseases. Diseases that cause structural disorders of the blood and affect the coagulation mechanism, such as leukemia and hemophilia.
- Inadequate bone volume and/or quality.
- Allergy or hypersensitivity to material.

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

**KULLANMA KILAVUZU /
USER MANUAL**

6. Warnings and Precautions

- Before the operation, the patient should be informed about the surgical risks and the positive and negative effects should be explained.
- The patient should be warned that the implanted device has a certain lifespan during which normal bone cannot replace, that it may be damaged by strenuous activity or trauma and may be replaced in the future.
- It is recommended to review the product-specific surgical technique before the surgeon performs the surgical procedure. IMPLATECH ONE IMPLANT can provide surgical technical information. Please contact the IMPLATECH ONE IMPLANT sales representative.
- Choosing the right implant is extremely important. Appropriate type, size, largest loss component is taken into consideration, the patient's age, activity levels, bone density are preferred depending on whether or not he has undergone any previous surgical operation.
- IMPLATECH ONE IMPLANT prostheses should be applied by securing respiration during oral use. Components or tools may injure the patient.
- Excessive bone loss or fracture of the dental implant may occur if an implant is loaded beyond its functional capacity. Physiological and anatomical conditions can affect the performance of dental implants.
- Misuse of small components inside the patient's mouth carries a risk of aspiration and/or swallowing.
- Forcing the implant deeper into the osteotomy than the depth created by the drills may result in implant, driver, or osteotomy damage.
- For short implants, clinicians should monitor patients closely for any of the following conditions: bone loss around the implant, changes in the response of the implant to percussion, or radiographic changes related to bone-to-implant contact along the length of the implant. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinician chooses a short implant, he or she should consider a two-stage surgical approach, attaching the short implant to an additional implant, and placing the widest possible fixture. The clinician should also wait longer for osseointegration and avoid immediate loading.
- Reuse of IMPLATECH ONE IMPLANT products labeled for single use may result in product contamination, patient infection, and/or device not performing as intended.
- It is recommended NOT to place implants with a diameter of less than 4 mm in the posterior regions.
- Dental Implant system elements should only be used with their own system elements.
- Implants should be placed in a suitable diameter, sufficient number and in an axis compatible with the dental arch.
- Burs suitable for implant size and diameter should be used.
- Against the risk of contamination, as soon as the implant is taken out of the sterile package, it should be placed in the slot opened with the help of the necessary equipment without touching the surface.
- Reuse of products labeled as disposable may result in product contamination, patient infection, and/or failure of the device to perform as intended. Therefore, they are not suitable for second use after use.
- The product is disposable. Do not reuse.
- Check the expiry date of the product before use.
- Do not use the product if a foreign substance or contamination is seen in the package or on the product.
- Do not use the product that has fallen to the ground.
- Since there is a risk of infection, dispose of the product as medical waste within the framework of legal procedures after use. IMPLATECH ONE IMPLANT implant and prosthetic parts are available in a variety of configurations. Abbreviations are used to assist on each product label.
- Should not be opened.
- These devices should only be used by trained professionals. The surgical and corrective techniques required for the proper use of these devices are largely specialized and complex procedures. Improper technique may result in implant failure, loss of supporting bone, fracture of the restoration, inhalation, ingestion and/or ingestion of loose screws. The option of immediate functional load may be considered when the clinician decides that adequate primary stability has been achieved.
- The following factors should be considered when placing dental implants: bone quality, oral hygiene and medical conditions due to blood disorders or uncontrollable hormonal disturbances. The recovery period can be varied depending on the bone quality at the implantation site, the tissue response of the implanted device, and the surgeon's assessment of patient bone density during the surgical procedure. Appropriate occlusion should be applied to the implant restoration to avoid applying excessive force to the implant during the healing process.

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

KULLANMA KILAVUZU / USER MANUAL

7. Complications And Side Effects

Dental implant treatment outcomes can be affected by many variables. Possible complications and side effects listed below may be observed after implant treatment.

- Local pain
- Micro hemorrhages
- Swelling
- Bruising
- Local inflammations
- Gingival injuries
- Difficulty biting, chewing and talking
- Systemic or local infections (including Periimplantitis, periodontitis, gingivitis, fistulas)
- Paresthesia
- Chronic pain due to nerve damage
- Damage to the opposing or adjacent teeth
- Bone damage
- Bone resorption
- Sinus perforations
- Mandible fractures
- Trismus
- Hypersensitivity and/or allergic reactions
- Other toxicity reactions
- Aspiration or swallowing of parts during operation
- Implant fractures
- Implant failures
- Mandatory removal of implants
- Failueres/fractures of prosthetic parts
- Undesirable aesthetic results
- Longer recovery times than expected

8. Compatibility Information

- Do not use implants supplied from different metals and different manufacturers and accessories other than IMPLATECH ONE IMPLANT together.
- IMPLATECH ONE IMPLANT dental implant systems and surgical sets are compatible with each other.

9. Sterilization

NOTE: Sterilization materials that you are sure of and that have appropriate biological indicators should be used.

- The implant is in sterile packaging. The products should be used sterile. The box cover should not be opened in any way before use. The implant has been sterilized by Gamma Irradiation. DO NOT AGAIN STERILE PRODUCTS!!
- Set parts are not sterile.

10. Storage and Transport

- Products should be kept in their original packaging. In case of damage to the sterile tube, the product should not be used. Products should be stored at 25+-5C, protected from direct sunlight. See labels for specific storage and handling rules.

11. Procedure

- Preparation of the implant site should be performed with the burs in the IMPLATECH ONE IMPLANT surgical set with a processing speed of 400-900 rpm under copious irrigation. Care should be taken not to overheat the bone. The use of burs that

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approverd By
ŞÜKRİYE AKARLAR

KULLANMA KILAVUZU / USER MANUAL

have lost their sharpness should be avoided. It can cause unnecessary trauma to the bone. For the drilling procedure, the Surgical Technique Guide should be consulted.

- Check the length and diameter of the implant on the product label.
- There is a label on the product packaging to ensure traceability. Affix this label to the patient's file. Open the outer pack and empty the sterile inner pack into the surgical field.
- Use the IMPLATECH ONE IMPLANT implant carrier to move the implant into the surgically prepared implant socket.
- IMPLATECH ONE IMPLANT abutments are placed on the relevant implants with the help of an internal screw and a hand wrench. A torque of 30-35 Ncm is recommended for optimal connection between parts.
- IMPLATECH ONE IMPLANT dental implant systems and surgical sets are compatible with each other.

12. Healing Process of Implants

- In general, implants should be left to heal for 2-4 hours, depending on the patient's bone quality, type, and general health. Immediate and early loading can be achieved by choosing the right case and adhering to the immediate loading protocols accepted in the literature. If immediate loading is to be done, the initial torque of the implant should not be less than 35Ncm, the occlusion environment should be controlled, if there is more than one implant, it should be splinted, and in cases of complete edentulism, at least 4 implants in the mandible and six implants in the maxilla should be connected to each other with prosthesis.

13. MRI Safety Information

IMPLATECH ONE IMPLANT Dental Implants have been found to be MR Conditional by testing in a non-clinical setting. Patients using this device can be safely scanned in an MR system when the following conditions are met:

- 1.5 T and 3.0 T static magnetic field
- Maximum surface area gradient of 3,000 gauss/cm (30 T/m)
- At maximum MR system, average specific absorption rate (SAR) of the whole body 2 W/kg (Normal Operating Mode)

Within the scanning conditions described above, IMPLATECH ONE IMPLANT Dental Implants are expected to have a maximum temperature rise of 4°C at 3.0 T and 3°C at 1.5 T after 15 minutes of continuous scanning.

In non-clinical testing, device-induced image artifact expands radially up to 2.7 cm and 2.2 cm from the implant when scanned with a gradient echo pulse array at 3.0 T and 1.5 T, respectively.

14. Disposal

Disposal of the products should be done in accordance with the rules within the framework of legal procedures and environmental requirements. Contaminated products and sharp parts are hazardous and should be disposed of as medical waste in appropriate conditions.

15. Information That Should Be Provided To The Patient

Patients should be informed about the contraindications, warnings, precautions, side effects and complications of Implatech One dental implants.

Patients should also be informed about MR safety information, and the diameter-height and serial number information of the implants used should be provided to the patient.

16. Disclaimer Of Liability

These products are part of a general concept and should only be used together with the relevant original products, in accordance with İmplatek Sağlık's instructions and recommendations. The use of products that are not distributed by İmplatek Sağlık and made by third parties will void any expressed or implied warranty or other obligations of İmplatek Sağlık. İmplatek Sağlık products should be used in accordance with the instructions for use provided by the manufacturer. It is the clinician's responsibility to use products in accordance with these instructions and to determine whether the product is appropriate for the individual patient situation. The clinician is also responsible for regularly following the latest updates on the relevant İmplatek Sağlık product and its applications. İmplatek Sağlık is not responsible for any direct, indirect, punitive or other damages arising out of or in connection with any errors in professional judgment or application related to the use of İmplatek Sağlık products; disclaims any liability, express or implied.

17. Additional Information

The shelf life of sterile products is 5 years.

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

**KULLANMA KILAVUZU /
USER MANUAL**

18. Symbols

Symbols and Meanings Prepared According to EN 15223-1: 2012 Standard

	Notified Body Number		Sterilized by Gamma Irradiation
	Do Not Sterilize For The Second Time		Lot number
	Do Not Use For The Second Time		Production date
	See User Manual		Expiration date
	Do Not Use If The Package Is Damaged		Do Not Expose To Direct Sunlight
	Manufacturer Information		Keep Away From Contact With Water
	Referance number		Barcode
	MRI Safety Information		Attention
	MR Koşullu		Dikkat, birlikte verilen belgelere danışın

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

KULLANMA KILAVUZU /
USER MANUAL

IMPLATECH ONE-S SHORT IMPLANT

**1. Product Description**

Implatech One-S Short Dental Implants are manufactured from biocompatible titanium or titanium alloy. Implatech One-S Short Dental Implants contain various surface treatments. Please refer to individual product labels for specific product description.

Implatech One-S Short dental implants are part of Implatech One-S Short Dental Implant Systems, which is an integrated system with related abutments, healing abutments, cover screws, surgical and prosthetic parts and instruments; They are endosseous screw type implants that have an internal conical octagon connection between the implant and the abutment. Surface roughening of Implatech One-S Short implants produced from biocompatible Titanium Grade 4 (ASTM F 67) material is done with biphasic calcium phosphate (BCP).

They are packaged together with the cover screw produced using Ti6Al4V-ELI (ASTM F 136) material.

Implatech One-S Short implants are produced in three different platforms as narrow platform (NP), regular platform (RP) Diameter and length options of Implatech One-S Short implants according to their platforms are given in the table below.

Platform	NP	RP	RP
Implant Ø(D)	Ø3.7	Ø4.1	Ø4.7
Uzunluk (L)			
6 mm	01.02.37.06	01.02.41.06	01.02.47.06

2. Intended Use

IMPLATECH ONE-S SHORT IMPLANT dental implants;

The implant, which is the main part of the dental implant system, is known as the artificial tooth root. It adheres effectively to the jawbone with its sandblasted and hydrophilic outer surface. Parts of Dental Implant Systems can be used to support and/or retain single-member, multi-member or all-arch fixed and/or removable dentures in partially or fully edentulous maxillary or mandibular arches. Implatech One-S Short implants are intended to be used in oral implantation as a support for prosthetic applications to be made in order to eliminate the function, phonation and aesthetic problems caused by tooth deficiencies in the lower and/or upper jaw bones.

3. Target Patient Group And Intended User

Implatech One-S Short implants are intended to be used in patients with complete or partial edentulism who have completed their growth and development and do not have the conditions specified in the contraindications.

Implatech One-S Short implants can only be used in atrophic jaws with inadequate bone height.

Implatech One dental implants are intended to be used by dentists only.

Clinicians must have sufficient knowledge of implantology and practice skills, as well as comply with the instructions for use, in order to use Implatech One dental implants safely and properly.

4. Indications

Implatech One-S Short dental implants are indicated for the functional, phonetic and aesthetic rehabilitation of the upper and/or lower jaw in fully or partially edentulous patients.

There is a wide range of indications, from the rehabilitation of single or multiple tooth deficiencies with fixed prosthetic applications in anterior and/or posterior regions where sufficient bone volume and quality is available, to the treatment of edentulous patients with fixed or removable prosthesis applications supported by implants.

After tooth extraction or loss, they can be applied as single or double-stage surgery with immediate, early or late implantation techniques.

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

**KULLANMA KILAVUZU /
USER MANUAL**

Implatech One-S Short implants are indicated for use in severely atrophic jaws where anatomical structures do not allow the placement of standard length implants due to insufficient bone height. It is indicated in cases where augmentations to increase the amount of vertical bone cannot be performed.

They are also indicated in terms of achieving primary stabilization with aggressive thread structure in atrophic bones with cancellous structure.

5. Contraindication

- Recent myocardial infarction. The first thing to consider here is surgical stress. Because the adrenaline that will be released during a surgical procedure will adversely affect the cardiac condition of the patient. In addition, there may be uncontrollable vasoconstriction, which may cause some disturbances in heart rhythm. Since almost all of these errors use anticoagulant drugs, a coagulation disorder may be in question. There is also a risk of infections.
- Those with a heart valve prosthesis. Likewise, there is a risk against surgical stress, coagulation imbalance and infections leading to the loss of the heart valve.
- Severe kidney diseases. Here, calcium is not absorbed from the tubules. This causes metabolic loss of calcium. Parathormone dysfunction. It causes metabolic osteopenia. risk of infection.
- Treatment of severe osteomalacia (rickets) There is hypophosphocalcic bone and osteoidosis. Failure to integrate the implant with the bone. 75% of rickets receive vit D + Ca. This hinders integration. There is a risk of infection.
- Generalized secondary osteoporosis. Although primary osteoporosis is physiological, secondary generalized osteoporosis has pathological conditions such as Hodking's disease. Rarefaction is seen in the bone structure. Absence of osteodosis; Although there is a mass increase in the bone, there is a volumetric decrease. Integration of the bone with the implant does not occur. There is a risk of infection.
- Uncontrolled Diabetes Mellitus. Hyperosmolarity is seen in the blood. Dehydration and metabolic diseases are in question. Angiopathy: Most diabetics suffer from micro and macro angiopathy. This leads to a predisposition to tissue degeneration. There is a risk of infection. In addition, wound healing is difficult.
- Those undergoing radiotherapy. The defense mechanism is broken. Osteoinduction and osteoconduction are impaired. Physiological periosteal activity is impaired. There is a predisposition to tissue necrosis and a risk of infection.
- Chronic or severe alcoholism. Liver diseases such as cirrhosis are in question. These, in turn, eventually lead to coagulation disorders. Medullary disorders (vit B1, B6, B12). As a result of these, there may be risk of anemia, thrombocyte disorders, hemorrhage in different distant places. Delay in recovery occurs due to malnutrition. Psychological disorders. There is a risk of infection.
- Severe hormonal disorders. There is a metabolic calcium deficiency. There is deterioration of the implant bed.
- Drug addiction. Loss of sense of pruritis Nutritional disorder. Loss of resistance to diseases. Psychological disorders. There is a risk of infection.
- Long-term use of immunosuppressant drugs. Delay in recovery. Medullary aplasia. Increased bone fragility. risk of infection
- Severe connective tissue diseases. For example, those that cause permanent damage to connective tissue, such as Lupus Erythematosus.
- Severe blood diseases. Diseases that cause structural disorders of the blood and affect the coagulation mechanism, such as leukemia and hemophilia.
- Inadequate bone volume and/or quality.
- Allergy or hypersensitivity to material.

6. Warnings and Precautions

- Before the operation, the patient should be informed about the surgical risks and the positive and negative effects should be explained.
- The patient should be warned that the implanted device has a certain lifespan during which normal bone cannot replace, that it may be damaged by strenuous activity or trauma and may be replaced in the future.
- It is recommended to review the product-specific surgical technique before the surgeon performs the surgical procedure. IMPLATECH ONE-S SHORT IMPLANT can provide surgical technical information. Please contact the IMPLATECH ONE-S SHORT IMPLANT sales representative.
- Choosing the right implant is extremely important. Appropriate type, size, largest loss component is taken into consideration, the patient's age, activity levels, bone density are preferred depending on whether or not he has undergone any previous surgical operation.

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

**KULLANMA KILAVUZU /
USER MANUAL**

- IMPLATECH ONE-S SHORT IMPLANT prostheses should be applied by securing respiration during oral use. Components or tools may injure the patient.
- Excessive bone loss or fracture of the dental implant may occur if an implant is loaded beyond its functional capacity. Physiological and anatomical conditions can affect the performance of dental implants.
- Misuse of small components inside the patient's mouth carries a risk of aspiration and/or swallowing.
- Forcing the implant deeper into the osteotomy than the depth created by the drills may result in implant, driver, or osteotomy damage.
- For short implants, clinicians should monitor patients closely for any of the following conditions: bone loss around the implant, changes in the response of the implant to percussion, or radiographic changes related to bone-to-implant contact along the length of the implant. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinician chooses a short implant, he or she should consider a two-stage surgical approach, attaching the short implant to an additional implant, and placing the widest possible fixture. The clinician should also wait longer for osseointegration and avoid immediate loading.
- Reuse of IMPLATECH ONE-S SHORT IMPLANT products labeled for single use may result in product contamination, patient infection, and/or device not performing as intended.
- It is recommended NOT to place implants with a diameter of less than 4 mm in the posterior regions.
- Dental Implant system elements should only be used with their own system elements.
- Implants should be placed in a suitable diameter, sufficient number and in an axis compatible with the dental arch.
- Burs suitable for implant size and diameter should be used.
- Against the risk of contamination, as soon as the implant is taken out of the sterile package, it should be placed in the slot opened with the help of the necessary equipment without touching the surface.
- Reuse of products labeled as disposable may result in product contamination, patient infection, and/or failure of the device to perform as intended. Therefore, they are not suitable for second use after use.
- The product is disposable. Do not reuse.
- Check the expiry date of the product before use.
- Do not use the product if a foreign substance or contamination is seen in the package or on the product.
- Do not use the product that has fallen to the ground.
- Since there is a risk of infection, dispose of the product as medical waste within the framework of legal procedures after use. IMPLATECH ONE-S SHORT IMPLANT implant and prosthetic parts are available in a variety of configurations. Abbreviations are used to assist on each product label.
- Should not be opened.
- These devices should only be used by trained professionals. The surgical and corrective techniques required for the proper use of these devices are largely specialized and complex procedures. Improper technique may result in implant failure, loss of supporting bone, fracture of the restoration, inhalation, ingestion and/or ingestion of loose screws. The option of immediate functional load may be considered when the clinician decides that adequate primary stability has been achieved.
- The following factors should be considered when placing dental implants: bone quality, oral hygiene and medical conditions due to blood disorders or uncontrollable hormonal disturbances. The recovery period can be varied depending on the bone quality at the implantation site, the tissue response of the implanted device, and the surgeon's assessment of patient bone density during the surgical procedure. Appropriate occlusion should be applied to the implant restoration to avoid applying excessive force to the implant during the healing process.

7. Complications And Side Effects

Dental implant treatment outcomes can be affected by many variables. Possible complications and side effects listed below may be observed after implant treatment.

- Local pain
- Micro hemorrhages
- Swelling
- Bruising
- Local inflammations
- Gingival injuries
- Difficulty biting, chewing and talking
- Systemic or local infections (including Periimplantitis, periodontitis, gingivitis, fistulas)
- Paresthesia
- Chronic pain due to nerve damage

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

KULLANMA KILAVUZU / USER MANUAL

- Damage to the opposing or adjacent teeth
- Bone damage
- Bone resorption
- Sinus perforations
- Mandible fractures
- Trismus
- Hypersensitivity and/or allergic reactions
- Other toxicity reactions
- Aspiration or swallowing of parts during operation
- Implant fractures
- Implant failures
- Mandatory removal of implants
- Failures/fractures of prosthetic parts
- Undesirable aesthetic results
- Longer recovery times than expected

8. Compatibility Information

- Do not use implants supplied from different metals and different manufacturers and accessories other than IMPLATECH ONE-S SHORT IMPLANT together.
- IMPLATECH ONE-S SHORT IMPLANT dental implant systems and surgical sets are compatible with each other.

9. Sterilization

NOTE: Sterilization materials that you are sure of and that have appropriate biological indicators should be used.

- The implant is in sterile packaging. The products should be used sterile. The box cover should not be opened in any way before use. The implant has been sterilized by Gamma Irradiation. DO NOT AGAIN STERILE PRODUCTS!!
- Set parts are not sterile.

10. Storage and Transport

- Products should be kept in their original packaging. In case of damage to the sterile tube, the product should not be used. Products should be stored at 25+-5C, protected from direct sunlight. See labels for specific storage and handling rules.

11. Procedure

- Preparation of the implant site should be performed with the burs in the IMPLATECH ONE-S SHORT IMPLANT surgical set with a processing speed of 400-900 rpm under copious irrigation. Care should be taken not to overheat the bone. The use of burs that have lost their sharpness should be avoided. It can cause unnecessary trauma to the bone. For the drilling procedure, the Surgical Technique Guide should be consulted.
- Check the length and diameter of the implant on the product label.
- There is a label on the product packaging to ensure traceability. Affix this label to the patient's file. Open the outer pack and empty the sterile inner pack into the surgical field.
- Use the IMPLATECH ONE-S SHORT IMPLANT implant carrier to move the implant into the surgically prepared implant socket.
- IMPLATECH ONE-S SHORT IMPLANT abutments are placed on the relevant implants with the help of an internal screw and a hand wrench. A torque of 30-35 Ncm is recommended for optimal connection between parts.
- IMPLATECH ONE-S SHORT IMPLANT dental implant systems and surgical sets are compatible with each other.

12. Healing Process of Implants

- In general, implants should be left to heal for 2-4 hours, depending on the patient's bone quality, type, and general health. Immediate and early loading can be achieved by choosing the right case and adhering to the immediate loading protocols accepted in the literature. If immediate loading is to be done, the initial torque of the implant should not be less than 35Ncm, the occlusion environment should be controlled, if there is more than one implant, it should be splinted, and in cases of complete

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

	93/42/EEC Tıbbi Cihaz Direktifi CE Teknik Dosya 93/42/EEC Medical Device Directive CE Technical File Product Name : Dental Implant System Substructures Ürün Adı : Dental İmplant Sistemi (Alt Yapılar)	Document No: TD.01/2.4.2 Release Date:02.01.2020 Revision No: 01 Revision Date: 24.04.2024 Sayfa 30 / 37
	KULLANMA KILAVUZU / USER MANUAL	

edentulism, at least 4 implants in the mandible and six implants in the maxilla should be connected to each other with prosthesis.

13. MRI Safety Information

IMPLATECH ONE-S SHORT IMPLANT Dental Implants have been found to be MR Conditional by testing in a non-clinical setting. Patients using this device can be safely scanned in an MR system when the following conditions are met:

- 1.5 T and 3.0 T static magnetic field
- Maximum surface area gradient of 3,000 gauss/cm (30 T/m)
- At maximum MR system, average specific absorption rate (SAR) of the whole body 2 W/kg (Normal Operating Mode)

Within the scanning conditions described above, IMPLATECH ONE IMPLANT Dental Implants are expected to have a maximum temperature rise of 4°C at 3.0 T and 3°C at 1.5 T after 15 minutes of continuous scanning.

In non-clinical testing, device-induced image artifact expands radially up to 2.7 cm and 2.2 cm from the implant when scanned with a gradient echo pulse array at 3.0 T and 1.5 T, respectively.

14. Disposal

Disposal of the products should be done in accordance with the rules within the framework of legal procedures and environmental requirements. Contaminated products and sharp parts are hazardous and should be disposed of as medical waste in appropriate conditions.

15. Information That Should Be Provided To The Patient

Patients should be informed about the contraindications, warnings, precautions, side effects and complications of Implatech One dental implants.

Patients should also be informed about MR safety information, and the diameter-height and serial number information of the implants used should be provided to the patient.

16. Disclaimer Of Liability

These products are part of a general concept and should only be used together with the relevant original products, in accordance with İmplatek Sağlık's instructions and recommendations. The use of products that are not distributed by İmplatek Sağlık and made by third parties will void any expressed or implied warranty or other obligations of İmplatek Sağlık. İmplatek Sağlık products should be used in accordance with the instructions for use provided by the manufacturer. It is the clinician's responsibility to use products in accordance with these instructions and to determine whether the product is appropriate for the individual patient situation. The clinician is also responsible for regularly following the latest updates on the relevant İmplatek Sağlık product and its applications. İmplatek Sağlık is not responsible for any direct, indirect, punitive or other damages arising out of or in connection with any errors in professional judgment or application related to the use of İmplatek Sağlık products; disclaims any liability, express or implied.

17. Additional Information

The shelf life of sterile products is 5 years.

Hazırlayan/ Prepared By ESEN ERDOĞAN 	Onaylayan/ Approved By ŞÜKRİYE AKARLAR 
--	--

**KULLANMA KILAVUZU /
USER MANUAL**

18. Symbols

Symbols and Meanings Prepared According to EN 15223-1: 2012 Standard

	2292	Notified Body Number		Sterilized by Gamma Irradiation
		Do Not Sterilize For The Second Time		Lot number
		Do Not Use For The Second Time		Production date
		See User Manual		Expiration date
		Do Not Use If The Package Is Damaged		Do Not Expose To Direct Sunlight
		Manufacturer Information		Keep Away From Contact With Water
		Referance number		Barcode
		MRI Safety Information		Attention

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Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approverd By
ŞÜKRİYE AKARLAR

IMPLATECH ONE IMPLANT COVER SCREWS

**NP Cover
Screws**

**RP Cover
Screws**

1. Product Description

Implatech One cover screws are part of Implatech One Dental Implant Systems, an integrated system with Implatech One dental implants and related abutments, healing abutments, surgical and prosthetic parts and instruments; It is produced as a single piece using Ti6Al4V-ELI (ASTM F 136) material.

Implatech One cover screws are components that cover the implant-abutment connection area during the healing process and prevent soft and hard tissue growth in this area.

While the upper part of the cover screw tightly covers the implant, the threaded part fits into the internal screw thread of the implant.

Implatech One cover screws are packaged with Implatech One dental implants and are supplied sterile.

Implatech One cover screws are compatible with all Implatech One dental implants.

Implatech One cover screws are designed in two different platforms, narrow platform (NP), normal platform (RP):

İmplant / Platform	Dar Platform (NP)	Normal Platform (RP)
Implatech One İmplant	Kapatma Vidası NP REF: 03.05.10.01	Kapatma Vidası RP REF: 03.05.20.01
Implatech One-S Short İmplant		

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

**KULLANMA KILAVUZU /
USER MANUAL**

2. Intended use

They are intended to used with the relevant Implatech One dental implants in the lower and/or upper jaws in partially or completely edentulous patients by covering the implant-abutment connection for preventing soft and hard tissue growth in this area during the healing process area in cases where two-stage surgical technique is used.

3.Target Patient Group And Intended User

Implatech One cover screws are intended for use with Implatech One dental implants in patients with complete or partial edentulism who have completed growth and development and do not have the conditions specified in the contraindications.

Implatech One cover screws are intended for use by dentists only.

Clinicians must have sufficient knowledge of implantology and practice skills, as well as comply with the instructions for use, in order to use Implatech One dental implants safely and properly.

4. Indications

Implatech One cover screws are indicated for use with Implatech One dental implants placed for the functional, phonetic and aesthetic rehabilitation of the upper and/or lower jaw in totally or partially edentulous patients.

5. Contraindications

The use of Implatech One cover screws is contraindicated in the presence of the following conditions;

- Patients who are medically unfit for implant treatment
- Allergy or hypersensitivity to Ti6Al4V-ELI (ASTM F 136) materials

6. Warnings / Cautions / Precautions

General

- Implatech One dental implant systems products should be applied by paying attention to the instructions for use provided by Implatek. It is the clinician's responsibility to use the products according to these instructions and to determine whether the implants are appropriate for the individual patient situation.
- Treatments with implants may lead to biological and mechanical problems, including bone loss and fatigue fracture of implants.
- Close collaboration between the surgeon, the prosthodontist and the dental laboratory technician is essential for successful implant treatment.
- Implatech One cover screws should only be used with their own system elements and paired with the correct platform during the surgical and prosthetic phases. The use of products with different brands and materials may lead to mechanical problems, failure of implants, tissue damage or aesthetic dissatisfaction.
- Implatech One cover screws are in sterile packaging. Do not resterilize or process the cover screws. Cleaning and sterilization can damage material and design features, leading to treatment failure.

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

**KULLANMA KILAVUZU /
USER MANUAL**

- Sterile handling is essential. Never use potentially contaminated components. Contamination can cause infections.
- All surgical instruments and tooling used must be kept in good condition and care must be taken not to damage the implants or other components of the instruments.

Pre-operative

- A comprehensive clinical and radiological examination is necessary to determine the psychological and physical condition of the patient.
- Preoperative hard and soft tissue deficiencies may result in undesirable aesthetic results and/or negative implant angulation.
- Particular attention should be paid to patients with local or systemic factors that may affect hard tissue healing, soft tissue healing, or osseointegration of implanted implants. (eg Type 1 Diabetes, bone metabolism diseases, bleeding disorders, anticoagulant therapy, parafunctional habits such as bruxism, smoking, poor oral hygiene, pregnancy, acute periodontal diseases, orofacial radiotherapy, adjacent tissue infections, etc.)
- Particular attention should be paid to patients receiving bisphosphonate treatment.
- Implant treatment should not be applied to patients who have not completed bone growth and development.
- Placement of implants and design of prostheses should be specific to the patient. In the presence of extreme bruxism and other parafunctional habits and jaw relationship disorders, different treatment options can be planned.

During the operation

- Before placing the Implatech One cover screws, ensure that the inner surface of the implant is clean and free from blood.
- Be careful against the risk of aspiration and swallowing due to the size of the products during intraoral application. Aspiration of products can lead to infection or unwanted physical injury.
- After implant placement, the clinician decides whether to use immediate or delayed loading protocols, based on the evaluation of the bone quality and primary stability in the operation.
- Quality and/or quantity of remaining bone, disturbed initial healing, local infections and general systemic diseases can be counted as causes of implant failures that occur immediately after the operation or after initial osseointegration has been achieved.
- Against the risk of contamination, it should be placed on the implant with the help of necessary equipment as soon as it is removed from the sterile package. Impairment of sterility may adversely affect the treatment outcome.

Post-operative

- For the long-term success of implant treatments, it is important to follow up the patient regularly and to give the necessary oral hygiene training to the patient.

7. Complications and side effects

Dental implant treatment outcomes can be affected by many variables. Possible complications and side effects listed below may be observed after the use of İmplatek cover screws.

- Hypersensitivity and/or allergic reactions
- Other toxicity reactions
- Aspiration or swallowing of parts during operation

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

KULLANMA KILAVUZU / USER MANUAL

- Triggering of the pharyngeal reflex during the insertion of the product
- Growth of bone tissue on the cover screws during the healing process
- In some cases, it becomes visible early as exposure during the healing process.

8. Compatibility information

Implatech One cover screws are compatible with İmplatek Sağlık implant systems products and components.

Do not use implants and components made of different metals and supplied from different manufacturers with İmplatek Sağlık dental implant system products.

9. Sterilization

Implatech One cover screws are delivered sterile and for single use only with Implatech One dental implants. They should not be cleaned or sterilized. İmplatek Sağlık does not accept any responsibility for the processing of devices that are initially delivered as sterile.

The protected packaging system protects the sterilized device from external factors and if the device is stored correctly, it maintains its sterility until the expiry date. When removing the device from its packaging, the rules of asepsis should be followed. The sterile package should not be opened until just before the insertion of the cover screw. Inspect the devices for damage before opening the sterile packaging. Products with damaged packaging system should not be used. A previously used or non-sterile cover screws should not be used under any circumstances.

10. Storage and Transport

The products should be stored at room temperature, in a dry place protected from direct sunlight. The products must be kept in their original packaging. Improper storage conditions may affect the product's characteristics and cause failure.

The expiry date of the product must be checked before using it. In case of damage to the sterile package, the product should not be used.

See labelings for specific storage and handling rules.

11. Procedure

- Take the Implatech One cover screw from the sterile package with the screwdriver.
- Connect the cover screw to the implant you placed in the socket and tighten it by hand.
- Loosen the Implatech One cover screw by hand using the screwdriver to remove it.

NOTE: When using the Implatech One cover screws, it is recommended to ensure that the screwdriver is firmly seated on the screw to avoid the risk of aspiration/swallowing of the screw.

NOTE: The recommended tightening torque is 5-10 Ncm manually by hand.

12. MRI Safety Information

Implatech One cover screws have been shown to be MRI conditional through testing in non-clinical tests. Patients using this device can be safely screened in an MRI system if the following conditions are met:

- 1.5 T and 3.0 T static magnetic field
- Maximum surface area gradient of 3,000 gauss/cm (30 T/m)
- At maximum MRI system, average specific absorption rate (SAR) of the whole body 2 W/kg (Normal Operating Mode)

Hazırlayan/ Prepared By
ESEN ERDOĞAN

Onaylayan/ Approved By
ŞÜKRİYE AKARLAR

	93/42/EEC Tibbi Cihaz Direktifi CE Teknik Dosya 93/42/EEC Medical Device Directive CE Technical File Product Name : Dental Implant System Substructures Ürün Adı : Dental İmplant Sistemi (Alt Yapılar)	Document No: TD.01/2.4.2 Release Date:02.01.2020 Revision No: 01 Revision Date: 24.04.2024 Sayfa 36 / 37
	KULLANMA KILAVUZU / USER MANUAL	

Under the scanning conditions described above, Implatech One dental implants are expected to have a maximum temperature rise of 4°C at 3.0 T and 3°C at 1.5 T after 15 minutes of continuous scanning.

In non-clinical tests, device-induced image artifact expands radially up to 2.7 cm and 2.2 cm from the implant when scanned with a gradient echo pulse array at 3.0 T and 1.5 T, respectively.

13. Disposal

Disposal of the products should be done in accordance with the rules within the framework of legal procedures and environmental requirements. Contaminated products and sharp parts are hazardous and should be disposed of as medical waste in appropriate conditions.

14. Information that should be provided to the patient

Patients should be informed about the contraindications, warnings, precautions, side effects and complications of Implatech One cover screws.

Patients should also be informed about MR safety information.

15. Disclaimer of liability

These products are part of a general concept and should only be used together with the relevant original products, in accordance with İmplatek Sağlık's instructions and recommendations. The use of products that are not distributed by İmplatek Sağlık and made by third parties will void any expressed or implied warranty or other obligations of İmplatek Sağlık.

İmplatek Sağlık products should be used in accordance with the instructions for use provided by the manufacturer. It is the clinician's responsibility to use products in accordance with these instructions and to determine whether the product is appropriate for the individual patient situation. The clinician is also responsible for regularly following the latest updates on the relevant İmplatek Sağlık product and its applications. İmplatek Sağlık is not responsible for any direct, indirect, punitive or other damages arising out of or in connection with any errors in professional judgment or application related to the use of İmplatek Sağlık products; disclaims any liability, express or implied.

16. Additional Information

The shelf life of sterile products is 5 years.

Hazırlayan/ Prepared By ESEN ERDOĞAN 	Onaylayan/ Approved By ŞÜKRİYE AKARLAR 
--	--

**KULLANMA KILAVUZU /
USER MANUAL**

17. Symbols

Symbols and Meanings Prepared According to EN 15223-1: 2012 Standard

	Notified Body Number		Sterilized by Gamma Irradiation
	Do Not Sterilize For The Second Time		Lot number
	Do Not Use For The Second Time		Production date
	See User Manual		Expiration date
	Do Not Use If The Package Is Damaged		Do Not Expose To Direct Sunlight
	Manufacturer Information		Keep Away From Contact With Water
	Reference number		Barcode
	MRI Safety Information		Attention

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