

# HOW TO CHOOSE CEREC LIBRARY

ImplatechOne is compatible with **Astratech Osseospeed Ev Library**.

Please find below compatibility chart for ImplatechOne platform for **CEREC**

## DIGITAL IMPRESSION TRANSFERS

### SCAN POST IMPRESSION TRANSFER

ImplatechOne			Plastic Connection	DENTSPLY ASTRA TECH OSSEOSPEED EV	
	Implant Ø	Scan Post Code		Implant Ø	Scan Post Code
<b>NP</b>	Ø3,3 - Ø3,7	NP-S Scan Post	S	Ø3,6	AT EV 3,6 S
<b>RP</b>	Ø4,1 - Ø4,7	RP-L Scan Post	L	Ø4,2	AT EV 4,2 L



## DIGITAL CAD/CAM

RESTORATION SOLUTIONS

### CEREC TI-BASE ENGAGED DIGITAL ABUTMENT

ImplatechOne			Plastic Connection	DENTSPLY ASTRA TECH OSSEOSPEED EV	
	Implant Ø	TiBase CODE		Implant Ø	TiBase CODE
<b>NP</b>	Ø3,3 - Ø3,7	NP-S TiBase H0.7 NP-S TiBase H2.5	S	Ø3,6	AT EV 3,6 GH1 S
<b>RP</b>	Ø4,1 - Ø4,7	RP-L TiBase H0.7 RP-L TiBase H2.5	L	Ø4,2	AT EV 4,2 GH1 L

